

EVO
Australia

EVO ONE

Low Height Overbed Table Compact Laminate Top



Side Rails

Easy grip, aluminium side rails.

Easy Locking Castors

Accessible from both sides.

Features & Advantages

EVO ONE Low Height range overbed table was designed specifically for use with low beds, with a minimum height of 640mm.

The unique design allows the table top to be easily lifted without having to operate the control handle and will remain at the position required with no further adjustments. The table top can only be lowered by operating the control handle. No bending is required.

The modern C-Base is extremely strong and durable and will fit in well with any room decor.

The table top is made of High Pressure Compact Laminate. Designed with infection control in mind, the flat surface is smooth, non-absorbent and free from gaps making it easy to clean. The hard wearing surface is impact resistance, as well as its ability to resist fading from harsh hospital cleaning chemical use, make it an excellent choice for Hospital.



Technical Specifications

Table top dimensions	900 x 400 mm
Table top	High Pressure Laminate (HPL) Or Thermoform
Frame	Powder Coated Aluminium & Steel
Height range	640 - 880mm
C-Base profile	90mm
Maximum load	15 kg (Evenly Distributed)
Castors	4 x 60mm, Nylon Twin Castors, 2 x Lockable
Colours	See colour sheet