


**EVO**  
Australia

# EVO ONE

## Tilt Top Overbed Table



### Side Rails

To prevent objects falling from the table top.

### Easy Locking Castors

Accessible from both sides.

### Features & Advantages

The unique design allows the table top to be easily lifted without having to operate the control handle and will remain at the position required with no further adjustments. The table top can only be lowered by operating the control handle. No bending is required.

With a large tilting section and a small fixed section for items that need to be kept level

The modern c-base is extremely strong and durable and will fit in well with any room decor.

The large thermoform top has been designed with infection control in mind, the flat surface is smooth, non-absorbent and free from gaps making it easy to clean.



### Technical Specifications

Table top dimensions	Tilt Area 662x460mm, Fixed Area 250x460mm
Table top	High Pressure Laminate (HPL) Or Thermoform
Frame	Powder Coated Aluminium & Steel
Height range	760 - 1020mm
C-Base profile	90mm
Maximum load	15 kg (Evenly Distributed)
Castors	4 x 60mm, Nylon Twin Castors, 2 x Lockable
Colours	See colour sheet