

EVO
Australia

EVO TWO

Two Way Bedside Cabinet



Features & Advantages

EVO TWO has been designed to be very strong and durable with High Pressure Compact Laminate panels and aluminium corners. The EVO TWO have been ergonomically designed to enable easy access from both sides for the patient and nurse. Locks for the top drawer are optional and can be provided with master key if required.

With infection control in mind all surfaces on the EVO TWO are smooth and non-absorbent making it easy to clean. EVO TWO made of HPL for durability and impact resistance, as well as its ability to resist fading from harsh hospital cleaning chemical use, make it an excellent choice for healthcare furniture.

**Easy Locking
Castors**
Accessible from
both sides.



Technical Specifications

Structure	Aluminium & High Pressure Laminate (HPL)
Overall Height	750mm
Overall Width	480mm
Overall Depth	480mm
Castors	4 x 70mm, Nylon Twin Castors, 2 x Lockable
Colours	See colours sheet